

# EXERCISES – BIOLOGICAL SIGNALS

Exercise 7 - SS 2014 – Michel Kana



# What will we do today?

---

1. **The physiology behind EEG**
2. **Structure of the EEG Signal**
3. **EEG measurement with BIOPAC**
4. **Summary**


# The physiology behind EEG

- ❑ The nervous system is composed of neurons and glial cells
- ❑ Neurons are excitable cells that generate and carry electrical signals that are called action potentials
- ❑ The electrical activity spreading through the head and reaches the scalp
- ❑ The resulting voltage differences can be recorded as the electroencephalogram (EEG)
- ❑ EEG reflects the summation of synchronous activity of many neurons with similar spatial orientations
- ❑ EEG offers a good millisecond temporal resolution, however with limited spatial resolution



# Structure of the EEG Signal

- **Characteristics of the EEG signal**
  - ▣ Frequency components up to 300 Hz.
  - ▣ Amplitude in order to  $\mu\text{volts}$ .



# Structure of the EEG Signal


## ❑ Characteristics of the EEG signal

- ❑ Frequency components up to 300 Hz.
- ❑ Amplitude in order to  $\mu$ volts.

## ❑ Waves in the EEG signal

- ❑ **Delta wave:** high amplitude neural oscillations in the frequency range of 0.5 – 4 Hz. They are primarily associated with deep sleep.
- ❑ **Theta wave:** neural oscillations in the frequency range of 4 – 8 Hz. They have been associated with access to unconsciousness, creative inspiration and deep meditation.
- ❑ **Alpha wave:** neural oscillations in the frequency range of 8 – 13 Hz. They are associated to a relaxed awareness without any attention or concentration. They should be reduced or eliminated by opening the eyes, by hearing unfamiliar sounds, by anxiety, mental concentration or attention.
- ❑ **Beta wave:** neural oscillations in the frequency range of 13 – 30 Hz. They are associated with active thinking, active attention, focus on the outside world, or solving concrete problems.
- ❑ **Gamma wave:** neural oscillations in the frequency range of 30 – 90 Hz. They are suggested to be related to consciousness.




# Structure of the EEG Signal



## □ Parameters in the EEG signal

- **Standard deviation (STD):** measure of activity in the brain wave
- **Average Value (AVG):** measure of mean amplitude
- **Cycle Count (CC):** measure of central frequency

# EEG Measurement



## □ Operational configuration

- Multichannel recording with an electrodes cap using 10 to 20 Ag–AgCl disks.
  - For example, C3 and C4 can be used to record the right and left finger movement related signals
- 1-Channel recording.
- Differential or referential electrodes setup.

## □ Noise sources

- Muscle movement artifact (0 to 1000 Hz)
- Motion artifact from electrode movements (0 to 20 Hz)
- Power line interference (60 or 50 Hz)

# Exercise 1: EEG measurement with BIOPAC

## ❑ Biopac MP35 measurement system


- EEG is recorded using Biopac SS2L wires plugged in the third channel.
- Electrodes are attached as depicted, avoiding hair between scalp and electrodes.

## ❑ Biopac Student Lab PRO software

- The acquisition is set up at a sampling rate of 500 Hz.
- Analog Channel CH1 should have the preset *Electroencephalogram EEG (.5-100Hz)*
- Calculation Channel C1 should have the preset *EEG alpha (8-13Hz)*
- Calculation Channel C2 should have the preset *EEG beta (13-30Hz)*
- Calculation Channel C3 should have the preset *EEG delta (0.5-4Hz)*
- Calculation Channel C4 should have the preset *EEG theta (4-8Hz)*
- Calculation Channel C5 should have the preset *EEG gamma (30-90Hz)*

## ❑ EEG parameters calculation

- Estimate the STD, AVG and CC of each of the five brain waves





# Exercise 2: Relaxed with eyes open without blinking

## □ Procedure

- Subject is instrumented for EEG measurement with Biopac.
- The subject should be seated with legs fully relaxed and keep eyes open, staring at the computer screen without blinking during min 10 seconds.

## □ Evaluation

- Estimate the STD of the all brain waves

# Exercise 3: Relaxed with eyes closed

## ☐ Procedure

- ☐ Subject is instrumented for EEG measurement with Biopac.
- ☐ The subject should be seated with legs fully relaxed and eyes closed during 120 seconds.

## ☐ Evaluation

- ☐ Estimate the STD, AVG and CC of the alpha wave
- ☐ Verify the following assertions:
  - ☐ Females tend to have higher mean frequency of alpha waves
  - ☐ Alpha amplitude tend to be higher in outgoing subjects

# Exercise 4: Mental math with eyes closed

## ☐ Procedure

- ☐ Subject is instrumented for EEG measurement with Biopac.
- ☐ The subject should be seated with legs fully relaxed and eyes closed.
- ☐ After 60 seconds of baseline measurement at rest, the subject should mentally find the remainder of 12345 divided by 12 during another 60 seconds.

## ☐ Evaluation

- ☐ Estimate the STD, AVG and CC of the alpha and beta waves
  - ☐ By how much did the amplitude of alpha waves changed?
  - ☐ Are there changes in the beta waves?

# Exercise 5: Hyperventilation with eyes closed

## □ Procedure

- Subject is instrumented for EEG and PPG measurement with Biopac.
- The subject should be seated with legs fully relaxed and eyes closed.
- The subject should increase breathing rate to 60 cycles per minute while breathing deeply during 120 seconds.

## □ Evaluation

- Estimate the STD, AVG and CC of the alpha wave
  - It is expected that carbon dioxide levels fall, pH increases, blood pressure decreases, overall brain activity increases, alpha rhythms increase. Can you verify?

# Team Projects

## ❑ Project 1: Cardiovascular Signal Analyzer

- ❑ Digital filtering of a raw PPG signal
- ❑ Extraction of PH (pulse height) and PP (peak-to-peak) values from a filtered PPG signal
- ❑ MAP estimation using PH
- ❑ Fourier transform of PP intervals and estimation of HF and LF
- ❑ Implementation in Matlab, if possible with an interactive GUI
  - ❑ User should be able to import the raw signal import from a Biopac text export
  - ❑ User should be able to enter the sampling frequency, signal type (ECG or PPG or both) and channel numbers
  - ❑ User should be able to filter the raw signal
  - ❑ User should be able to execute PP, PH, MAP, LF, HF computation
  - ❑ User should be able to display plots of the raw signal for a given start and end timestamp
  - ❑ User should be able to display plots of PP, PH, MAP over the time for a given start and end timestamp and display the value of LF and HF

# Team Projects

## ❑ Project 2: Nervous Activity Analyzer

- ❑ Digital filtering of a raw EEG signal
- ❑ Extraction of alpha, beta, theta, delta waves from a filtered EEG signal
- ❑ Computation of STD, AVG and CC
- ❑ Implementation in Matlab, if possible with an interactive GUI
  - ❑ User should be able to import the raw signal import from a Biopac text export
  - ❑ User should be able to enter the sampling frequency
  - ❑ User should be able to filter the raw signal
  - ❑ User should be able to execute alpha, beta wave, theta, delta wave computation using Fourier or Wavelet transform or digital filtering
  - ❑ User should be able to execute STD, AVG, CC computation
  - ❑ User should be able to display plots of the raw signal for a given start and end timestamp
  - ❑ User should be able to display plots of alpha, beta wave, theta, delta waves over the time for a given start and end timestamp and display the values for STD, AVG and CC

# Team Projects

## ❑ **Project 3: Muscle Activity Analyzer**

- ❑ Digital filtering of a raw EMG signal
- ❑ Computation of rectified EMG from a filtered EMG signal
- ❑ Computation of the spectrum of the filtered EMG signal using Fourier transform
- ❑ Computation of RMS, ARV
- ❑ Implementation in Matlab, if possible with an interactive GUI
  - ❑ User should be able to import the raw signal import from a Biopac text export
  - ❑ User should be able to enter the sampling frequency
  - ❑ User should be able to filter the raw signal
  - ❑ User should be able to execute rectified EMG computation
  - ❑ User should be able to execute Fourier transform of the rectified EMG for a given start and end timestamp
  - ❑ User should be able to execute RMS, ARV computation for a given start and end timestamp
  - ❑ User should be able to display plots of the raw EMG, rectified EMG, EMG Fourier transform for a given start and end timestamp
  - ❑ User should be able to display the values for RMS, ARV for a given start and end timestamp