


EXERCISES – BIOLOGICAL SIGNALS

Exercise 1 - SS 2014 – Michel Kana

What will we do today?

1. **Semester plan**
2. **Pre-test**
3. **The physiology behind ECG**
4. **Structure of the ECG Signal**
5. **ECG measurement with BIOPAC**
6. **Heart rate processing with Matlab**
7. **Summary**

Semester plan



Pre-test


The physiology behind ECG

Electrical conduction in the heart


- The sinoatrial node (SA node) automatically depolarizes.
- Electrical activity goes rapidly to atrioventricular node (AV node) via internodal pathways.
- Depolarization spreads slowly across atria. Conduction slows through AV node.
- Depolarization moves rapidly through ventricular conducting system to the apex of the heart.
- Depolarization wave spreads upward from the apex.

The Electrocardiogram

- The ECG reflects the electrical activity of the heart.
- An ECG tracing shows the summed electrical potentials generated by all cells of the heart.
- The Einthoven's triangle defines electrodes position for 3-lead ECG measurement.



Ref.: D. Unglaub, Human Physiology



Structure of the ECG Signal

□ Waves in the ECG


- **P wave:** corresponds to the depolarization of atria. Duration 0.06 - 0.11 sec.
- **QRS complex** (Q, R and S waves): represents the progressive wave of ventricular depolarization. Duration 0.08 - 0.12 sec.
- **T wave:** represents the repolarization of ventricles.
- **U wave:** represents the repolarization of the papillary muscles or Purkinje fibers.

□ Intervals in the ECG

- **P-R interval:** begins at the onset of the P wave and ends at the onset of the QRS complex and represents the time the impulse takes to reach the ventricles from the sinus node. Duration 0.12 – 0.20 sec.
- **Q-T interval:** begins at the onset of the QRS complex and ends at the endpoint of the T wave and represents the duration from the depolarisation to the repolarisation of the ventricles. Duration 0.35 – 0.43 sec.
- **R-R interval:** time measurement between the R wave of one heartbeat and the R wave of the preceding heartbeat.

□ Segments in the ECG

- **P-R segment:** begins at the endpoint of the P wave and ends at the onset of the QRS complex and represents the duration of the conduction from the atrioventricular node, down the bundle of His and through the bundle branches to the muscle.
- **S-T segment:** begins at the endpoint of the S wave and ends at the onset of the T wave. During the ST segment, the atrial cells are relaxed and the ventricles are contracted so electrical activity is not visible.




Ref.: <http://noodle.med.yale.edu>

ECG Measurement

Electrodes and Connecting wires

- Electrolyte paste with chloride ions for better signal transmission
- Bipolar electrodes: record the voltage differential between the wrists and the legs.
- Unipolar electrodes: record the voltage difference between a reference electrode and the body surface.
- Plate electrodes: use metal disks held onto the skin with adhesive tape.
- Suction electrodes: use metal disks with a vacuum system to remain in place.
- Fluid column electrodes: designed to avoid direct contact with the skin
- Flexible electrodes: designed for infants with fine stainless steel or silver wire, attached to the skin like a small bandage.



Ref.: <http://www.edn.com>

Amplifier

- Converts the weak electrical signal from the body into a more readable signal.
- Includes optical isolators, preventing the possibility of accidental electric shock.

Storage and transmission device

- Includes built-in analog to digital converter.



Ref.: <http://www.biopac.com>

Biopac Student Lab PRO software

The screenshot displays the Biopac Student Lab PRO software interface. The window title is "Biopac Student Lab PRO® - [C:\Program Files\Biopac Student Lab\BSL 3.7\BSL PRO 3.7...". The menu bar includes File, Edit, Transform, Display, Window, MP35, and Help. The toolbar contains various icons for file operations, data manipulation, and display settings. Below the toolbar, there are measurement boxes showing "mean 7.38094 Kg", "max 23.96066 Kg", and "delta T 22.33200 sec". The channel boxes show "1 mean" and "3 max". The marker tools section includes a text input field with "Increasing grip strength - left hand, 11:01:52 AM" and a dropdown menu. The graph window displays two channels: "Force" (blue line) and "EMG" (red line). The Force channel has a vertical scale from -4.00 to 28.00 Kg, and the EMG channel has a vertical scale from -3.00 to 3.00 mV. The horizontal axis is labeled "seconds" and ranges from 0.00 to 112.00. The status bar at the bottom shows "Ready, Grid Off", "Chart", and "Disk".

File Name
Menu Bar
Toolbar

Channel Boxes
Append Markers
Marker Label
Event Markers

Channels

Horizontal (Time)
Scale and Scroll Bar

Journal Tools

Journal

Status Bar

Measurement Boxes
Channel | Measurement | Result

Marker Tools

Selected Area

Vertical (Amplitude)
Scales and Scroll Bar

Selection Tools
(Arrow, I-Beam, Zoom)

Status Indicator

Data, or Graph Window

Journal Window


ECG measurement with BIOPAC

□ Biopac MP35 measurement system

- 3 leads ECG is recorded using Biopac SS2L connecting wires plugged into the first channel.
- Bipolar plate electrodes with electrolyte paste are used.

□ Biopac Student Lab PRO software

- The Biopac MP35 system is connected to a PC where the Biopac Student Lab PRO software is installed for continuous signal acquisition
- The acquisition is set up for *Recording and Append into Memory* at a sampling rate of 500 Hz. The total acquisition length can be set to 2 hours.
- In the *Setup Channels* menu, the check boxes *Acquire Data*, *Plot on Screen*, *Enable Value Display* are enabled for the first channel.
- Analog Channel CH1 should have the preset ECG (.5 - 35 Hz). The frequency range of an ECG signal is of 0.05–250 Hz and its dynamic range of 1–10 mV
- Calculation Channel C1 should have the preset *R-R Interval*
- Calculation Channel C2 should have the preset *Heart Rate*



Exercise 1

❑ Sampling rate effect on the ECG

- ❑ MP menu > Set up Acquisition
- ❑ No significant statistical difference between 500 Hz and 250 Hz
- ❑ Undersampling (< 125 Hz) shows significant reductions in peak amplitude values
- ❑ 512 Hz ensures that the QRS complex is fully digitized. $512 \text{ Hz} = \text{sampling period of approximately } 1/512 = 2 \text{ ms}$, which captures at least $0.08 \text{ s} / 2 \text{ ms} = 40$ samples are captured during the shortest QRS complex.

❑ Hardware Settings for ECG

- ❑ MP menu > Set up Channels > View/Change Parameters
- ❑ Magnitude amplification from $1 - 2 \text{ mV}$ to about $0.5 - 1 \text{ V}$ with a total gain of about 500
- ❑ Add or subtract a constant to the signal in case the electrode has inherent offset
- ❑ AC Coupling is appropriate for ECG

Exercise 2

Hardware-based filters in Biopac

- MP menu > Set up Channels > View/Change Parameters
- Within the MP35/30 hardware itself before the data is displayed on your computer
- Cannot be used to filter data below 30 Hz
- Three programmable digital second-order filters: low pass, high pass, band pass, band stop.
- Two selectable hardware filter: a switchable 0.05 Hz / 0.5 Hz only for AC coupled signals, a switchable 1 kHz / 5 kHz low pass filter

Application: noise filtering in ECG

- Magnitude amplification from 1 – 2 mV to about 0.5 – 1 V with a total gain of about 500.
- High-pass filter (cut-off 0.03 to 5 Hz, usually 0.05 Hz) for the removal of low-frequency noise produced by respiration and electrode movement that results in a base line drift of the ECG signal
- Motionless subject to avoid EMG noise which has a wide frequency range (1 – 5000 Hz)
- Low-pass (cut-off 35 to 60 Hz) or notch filter for the removal of 60 Hz noise of power line interference
- Q setting: an optimally dampened filter has Q of 0.707 for low and high-pass filters, 5.000 for band pass and band stop filters

Input Channel Parameters

Channel Number: CH1

Channel Preset: CH1 Input

Channel Label: CH1 Input

Digital Filters:

Filter: 1

Type: Low Pass

Freq: 100.0000 Hz

Q: 0.707000

Hardware:

Gain: x200

Offset: 0. mV

Input coupling:

☐ AC ☒ DC

☐ 0.05Hz HP ☐ 0.5Hz HP

☐ 1 kHz LP ☒ 5 kHz LP

New Channel Preset Scaling... Cancel OK

Exercise 3

Calculation Channel Presets

- Rate Calculation: extracts information about the interval between a series of peaks in a waveform.
 - Function *Interval* for inter-beat interval calculation (i.e. RR interval durations)
 - Function *Rate (BPM)* for beats per minutes calculation
 - Function *Peak Maximum* for systolic blood pressure estimation from R-wave
 - Function *Peak Minimum* for diastolic blood pressure estimation from R-wave
 - Function *Mean Rate Value* for mean arterial pressure estimation from R-wave
- Filter Calculation: performs real-time digital filtering on Analog or Calculation channels.
 - Low Pass filter: pass only frequencies below the cutoff point
 - High Pass filter: pass only frequencies above the cutoff point
 - Band Pass (low + high) filter: pass only frequencies in the range specified by a low and high cutoff
 - Band Pass filter: pass only a narrow band of frequencies around a center cutoff
 - Band Stop filter: pass only a wide band of frequencies except of those around a center cutoff. E.g. for removing 60 Hz power line noise from ECG.
 - Band Stop – Line Freq filter: similar to Band Stop for 50 Hz or 60 Hz cutoff.

Off

New

ECG - R-R Interval

ECG - R-wave Amplitude

EEG alpha (8 - 13 Hz)

EEG beta (13 - 30 Hz)

EEG delta (1 - 5 Hz)

EEG theta (4 - 8 Hz)

EEG gamma (30 - 90 Hz)

EGG (.02 - .125 Hz)

EMG - Integrated (estimate)

EMG - RMS

Heart Rate (from ECG)

Lung Volume

Pulse Rate (from PPG)

Respiration Rate

dp/dt @ 200 samples/sec.

dp/dt @ 500 samples/sec.

dp/dt @ 1000 samples/sec.

Large Animal dp/dt Minimum

Large Animal dp/dt Maximum

Large Animal Systolic Blood Pressure (BP)

Large Animal Diastolic Blood Pressure (BP)

Large Animal Mean Blood Pressure (BP)

Large Animal Heart Rate (from ECG)

Small Animal dp/dt Minimum

Small Animal dp/dt Maximum

Small Animal Systolic Blood Pressure (BP)

Small Animal Diastolic Blood Pressure (BP)

Small Animal Mean Blood Pressure (BP)

Small Animal Heart Rate (from ECG)

Integrate

Smoothing

Difference

Rate

Math

Function

Filter

Expression

Delay

Exercise 4

□ Deep Breath Test

□ Materials

- A chronometer to measure inspiration and expiration duration.

□ Procedure

- The examiner explains to the subject how the test will be conducted: *“subject will be taught to breathe at a rate of 6 respiration cycles per minute, which corresponds to 0.1 Hz: 5 seconds for each inhalation and 5 seconds for each exhalation. The examiner will pace the breathing with a chronometer and give a vocal signal for the start of each inhalation and exhalation”*.
- Subject is instrumented by the examiner.
- Operator starts recording in Biopac Student Lab PRO software and observes heart rate. When both signals have stabilized for at least 60 seconds he informs the examiner to start with the test.
- Examiner gives a vocal signal to the subject to inhale and starts the chronometer.
- Operator adds a DEEP INHALATION event in Biopac Student Lab PRO software.
- Subject performs a 5 seconds inhalation.
- After the chronometer has run for 5 seconds, Examiner gives a vocal signal to the subject to exhale.
- Operator adds a DEEP EXHALATION event in Biopac Student Lab PRO software.
- Subject performs a 5 seconds exhalation.
- The previous steps are repeated 6 times.
- Operator continues recording in Biopac Student Lab PRO software and observes heart rate and respiration rate. When both signals have stabilized he stops recording after another 60 seconds.
- Operator exports the measurement File > Save As > *.TXT > Option.

Exercise 5

□ Deep Breath Test

▣ Data import from Biopac into Matlab

- Create a Matlab variables for the heart rate signal.

▣ Automated evaluation in Matlab

- Write a Matlab function for calculating the following values:
 - The *expiratory-to-inspiratory difference* ($E-I_{mean}$) is calculated as the mean of differences between the maximal heart rate during inspiration and minimal heart rate during expiration for all six breathing cycles.
 - The *expiratory-to-inspiratory ratio* (E/I_{mean}) is calculated as the mean of ratios of the maximal heart rate during inspiration divided by minimal heart rate during expiration for all six breathing cycles.

Summary

[Semester plan]

13 exercise sessions.

You can get 10 points for your attendance, 10 points from the first test and 10 points from the final presentation.

[ECG Physiology]

Electrical conduction in the heart and the origin of the electrocardiogram

[Structure of the ECG]

Waves, segments and waves in a normal electrocardiogram.

[ECG measurement and processing]

ECG recording with Biopac MP35. Hardware filters. Calculation and extraction of ECG parameters. Postprocessing in Matlab.

[Plan for the next week]

Photoplethysmogram, Galvanic skin response measurement & analysis.

Experiments with active change of posture. Bring materials with you!

