

EXERCISES – BIOLOGICAL SIGNALS


Exercise 4 - SS 2014 – Michel Kana

What will we do today?

1. **The physiology behind PPG**
2. **Structure of the PPG Signal**
3. **PPG measurement with BIOPAC**
4. **Summary**

The physiology behind PPG

- Photoplethysmography (PPG) is an optical measurement technique that can be used to detect blood volume changes in the microvascular bed of tissue
- A photodetector can be placed on a tissue's surface alongside the light-emitting diode and record the light that returns back
- Light intensity is attenuated by oxygenated and deoxygenated hemoglobin (in blood cells), myoglobin (in muscle), and cytochromes, as well as melanin (in skin)
- Reductions in light intensity indicate relative increases in blood volume and vice versa




Wiggers diagram

Structure of the PPG Signal


□ Waves in the PPG

- The photoplethysmogram waveform consists of one systolic wave and one diastolic wave.
- The second derivative photoplethysmogram waveform consists of four systolic waves (a, b, c, and d waves) and one diastolic wave (e wave).




Ref.: M. Elgendi "Standard Terminologies for Photoplethysmogram Signals"

Structure of the PPG Signal



Parameters in the PPG signal



- PH: pulse height or pulse amplitude
- CP: cardiac period
- FWHM: full width half max
- PW: peak width
- PT: peak threshold
- PPT: pulse travelling time





Information in the PPG

- Heart rate: $\frac{60}{CP}$ in beats per minute
- Blood pressure:
 - $p_1 - p_2 \cdot PTT$ (e.g. $p_1 = 246, p_2 = 0.4, PTT=315\text{ms}$)
 - $p_1 \cdot e^{p_2 \cdot PH} + p_3 \cdot e^{p_4 \cdot PH}$ (e.g. $p_1 = 105, p_2 = -4, p_3 = 0.2, p_4 = 18$)
- Respiration: using low frequency respiratory induced intensity variations and PPT
- Arterial stiffness: $\frac{c}{a}$
- Vascular aging: $\frac{b-c-d-e}{a}$

PPG Measurement



Ref.: <http://www.engin.swarthmore.edu/~dluong1/E72/FinalProject/heart.htm>



Operational configuration

- ☐ Side-by-side mode
- ☐ Transmission mode

Noise sources

- ☐ Movement artefacts
- ☐ Ambient light interference
- ☐ Variation in temperature
- ☐ power line interference

Exercise 1: PPG measurement with BIOPAC

❑ Biopac MP35 measurement system

- ❑ PPG is recorded using Biopac SS4LA plethysmograph transducer plugged in the first channel.
- ❑ Transducer should be attached to the left hand, index finger.

❑ Biopac Student Lab PRO software

- ❑ The acquisition is set up for *Recording* and *Append into Memory* at a sampling rate of 200 Hz. The total acquisition length can be set to 1 hour.
- ❑ In the *Setup Channels* menu, the check boxes *Acquire Data*, *Plot on Screen*, *Enable Value Display* are enabled for the first channel.
- ❑ Analog Channel CH1 should have the following settings
 - ❑ AC coupling 0.5 Hz High Pass, 1K Low Pass
 - ❑ Hardware-based filter 60 Hz Low Pass
 - ❑ Total gain 5000
- ❑ Calculation Channel C1 should calculate the Difference of the channel CH1 in order to estimate the first derivative of the PPG waveform. Channel C1 should not be visible.
- ❑ Calculation Channel C2 should calculate the Difference of the channel C1 in order to estimate the second derivative of the PPG waveform

❑ PPG parameters calculation

- ❑ Select two consecutive cardiac cycles
- ❑ Read the values for PH, CP, FWHM, PW, PT
- ❑ Estimate the heart rate
- ❑ Estimate arterial stiffness and vascular aging

Exercise 2: ECG & PPG measurement with BIOPAC

❑ Biopac MP35 measurement system

- ❑ PPG is recorded using Biopac SS4LA plethysmograph transducer plugged in the first channel.
- ❑ 3 leads ECG is recorded using Biopac SS2L connecting wires plugged in the 2nd channel.

❑ Biopac Student Lab PRO software

- ❑ Analog Channel CH2 should have the following settings
 - ❑ AC coupling 0.5 Hz High Pass, 1K Low Pass
 - ❑ Hardware-based filter 35 Hz Low Pass
 - ❑ Total gain 2000
- ❑ Calculation Channel C3 should calculate the heart rate from the PPG signal
- ❑ Calculation Channel C4 should calculate the heart rate from the ECG signal
- ❑ Calculation Channel C5 should estimate the mean arterial pressure from the pulse amplitude in the PPG signal

❑ Pulse travel time calculation

- ❑ Select one cardiac cycle
- ❑ Use the measurement window and calculate the PPT
- ❑ Estimate the mean arterial blood pressure using the PPT and the pulse amplitude
- ❑ Suggest a method for calculating and displaying the PPT during realtime measurement

❑ Heart rate comparison

- ❑ Select a sufficiently long time period
- ❑ Calculate and compare the mean, min and max heart rates obtained from ECG and PPG

Exercise 3: active standing test

□ **Active Change of Posture Test**

- Subject is instrumented for PPG measurement with Biopac by the examiner.
- The examiner records the PPG signal for 120 seconds.
- The examiner gives a vocal signal to the subject to stand.
- Subject stands up.
- The examiner records the PPG signal for 120 seconds.

■ **Evaluation**

- Identify significant phases in the heart rate and blood pressure waves
- Compare the mean values between phases

Summary

[What did we learn today]

Physiology behind the photoplethysmogram.

Structure of the PPG signals.

Measurements and signal analysis.

[Plan for the next week]

Respiration measurement & analysis