

Nervous Activity Analyzer

By Vicktoria Antonova & Anne Leijsen
Biological Signals, May 2013


Overview of Presentation

- Our Task
- Graphical User Interface (GUI)
- Matlab code
- Discussion
- Conclusion

Our Task:

- ❑ Digital filtering of a raw EEG signal
- ❑ Extraction of alpha, beta, theta, delta waves from a filtered EEG signal
- ❑ Computation of STD, AVG and CC
- ❑ Implementation in Matlab, if possible with an interactive GUI
 - ❑ User should be able to import the raw signal import from a Biopac text export
 - ❑ User should be able to enter the sampling frequency
 - ❑ User should be able to filter the raw signal
 - ❑ User should be able to execute alpha, beta wave, theta, delta wave computation using Fourier or Wavelet transform or digital filtering
 - ❑ User should be able to execute STD, AVG, CC computation
 - ❑ User should be able to display plots of the raw signal for a given start and end timestamp
 - ❑ User should be able to display plots of alpha, beta wave, theta, delta waves over the time for a given start and end timestamp and display the values for STD, AVG and CC

Graphical User Interface



Matlab Code

The screenshot displays the MATLAB development environment. The main window is the Editor, showing a script named `Browserfile.m`. The script contains two functions: `pushbutton1_Callback` and `column_value_Callback`. The `pushbutton1_Callback` function handles a button press by opening a file selector, reading the file, and plotting the data. The `column_value_Callback` function handles a column selection event. The script also includes comments and variable assignments for standard deviation and mean amplitude.

```
75  
76 % --- Executes on button press in pushbutton1.  
77 function pushbutton1_Callback(hObject, eventdata, handles)  
78 % hObject      handle to pushbutton1 (see GCBO)  
79 % eventdata    reserved - to be defined in a future version of MATLAB  
80 % handles      structure with handles and user data (see GUIDATA)  
81 [filename pathname] = uigetfile({'*.txt'}, 'File Selector');  
82 fullpathname = strcat(pathname, filename);  
83 text = fileread(fullpathname);  
84 set(handles.text2, 'String', fullpathname) % showing FullPathName  
85  
86 column = str2num(get(handles.column_value, 'string')); % get the column number and convert to num  
87  
88 Fs = str2num(get(handles.freq_value, 'string'));  
89 t1_value = str2num(get(handles.t1_value, 'string'));  
90 t2_value = str2num(get(handles.t2_value, 'string'));  
91 t1 = t1_value*Fs+1;  
92 t2 = t2_value*Fs;  
93  
94 x = importdata(fullpathname); % all data imported from file  
95 data = x(t1:1/Fs:t2,column); % take only data from column of interest  
96 plot(data);  
97  
98 s = std(data); % Find standard deviation  
99 set(handles.text12, 'String', s) % Write standard deviation  
100 mean_amplitude = mean(data); % Find mean amplitude  
101 set(handles.text13, 'String', mean_amplitude); % Write mean amplitude  
102  
103 function column_value_Callback(hObject, eventdata, handles)  
104 % hObject      handle to column_value (see GCBO)  
105 % eventdata    reserved - to be defined in a future version of MATLAB
```

The Workspace panel on the right shows the current workspace variables. The Command History panel shows the commands entered in the Command Window.

Name	Value
guide	
Browserfile	
guide	
Browserfile	
guide	
Browserfile	
guide	
Browserfile	
guide	
Browserfile	
help data	
Browserfile	

Discussion

- Program is made for .txt files (exported e.g. from Biopac), this should be taken into account using the program.
 - > Data should be properly formatted (use of .).
 - > Mention the right data column.
- Button refers to SPTool to change filter design (advanced users).
- Cycle Count is not showed, no wavelet transform.

Conclusion

- Program is working for available data.
- With future improvement, the Nervous Activity Analyzer could be used as part of e.g. ambulatory EEG recording.

Nervous Activity Analyzer

The new user interface for EEG recording

Thank you for your attention.
Are there any questions?